

Patented
Berkeley
Nitric Oxide
Test Strips




Berkeley Test®

A healthy
cardiovascular
system

How:
5 sec
1. It's simple,
place the
test strip with
the "saliva here" side
on your tongue.

A Nobel Prize was awarded for the discovery that Nitric Oxide (NO) plays a significant role in the cardiovascular system. NO is naturally made in your body, but as you age you produce less of it. Berkeley Life helps support NO levels for Heart Health by supplementing your diet with dietary nitrates. The Patented Berkeley Test Strips are a non-invasive saliva test. A diet rich in vegetables with high nitrate levels is best (see Berkeleylife.com for more information) or try The Berkeley Life Supplements. They provide the dietary nitrates found in plants, as well as key vitamins, to help support NO levels, cardiovascular health and overall health and vitality.

Directions:
You can test at any time during the day to check your NO levels. Ensure you have not consumed anything 10 minutes prior to testing. Allow a minimum of 1 to 2 hours after your nitrate high meal or supplement before testing.



Nitric Oxide (NO) Saliva Test Strips

- 10 sec**
2. Fold the strip over and gently press the two sides together.
 3. Now simply separate and compare color on test pad to the NO scale.
 4. Download the Berkeley Test app to get your exact reading and to track over time.



STRIPS ARE FOR SINGLE USE ONLY. DO NOT PUT TEST PAD END IN MOUTH.

15 sec
Self test
in seconds


WARNING:
KEEP OUT OF REACH OF CHILDREN. STORE IN A DRY PLACE AND AVOID EXCESSIVE HEAT.

Customer Care Toll Free:
844-512-0943 • berkeleylife.com
Manufactured for and distributed by:
L2G North America, Inc., 355 N. Canal
Street, Chicago, IL 60606, USA



10 Nitric Oxide Test Strips





Lifes2good

Proof No.: 02

Date: 14/10/2016

Part Code: PKG-BL-NOS-GL-CAR-1700.3

Artwork Title: LG2179.02_Berkeley Test Strips_Carton_10 Pack

Brand: Berkeley Life Artwork Item: Carton

UPC Barcode (80%): 817105020181 Country: USA

Size: 112mm x 40mm x 27mm

Printing Colours:

Cyan	Magenta	Yellow	Black	P 1505
P 705	P 707	P 708	P 212	P 214

Non-Printing Colours:

Keyline
