



keto mojo™

β-ketone
test strips

10

for use with keto-mojo β-ketone &
blood glucose monitoring system only.

what's your
mojo?



Apply blood to tip of strip only.
Sample will be pulled into strip.



Do not apply blood to strip surface.
Holding or touching strip will not
affect reading.

Chemical components in sensor:

- | | |
|--|---------|
| 1. β-Hydroxybutyrate Dehydrogenase (Pseudomonas sp.) | ≥ 0.5 U |
| 2. Mediator | 55% |
| 3. NAD | ≥ 0.5 g |
| 4. Enzyme protector | 8% |
| 5. Non-reactive ingredients | 29% |



Distributed by
KETO-CHECK, INC.
952 School St. Suite 212
Napa, CA 94559-2824
Tel: 800.513.1965
10am-4pm PST

for self-testing



302-4279301-013