

USAGE: As a dietary supplement, take 2 capsules per day with a meal or a glass of water or as directed by your qualified healthcare professional.

WARNING: Do NOT use this product if you are pregnant, lactating or trying to conceive. **NOTE:** If you have a medical condition, are under the age of 18, or taking medications, consult your healthcare professional before using this product.

Red Yeast Rice – Rice that has been fermented by the yeast *Monascus purpureus* has been used in Asia for centuries as a food preservative, food colorant and flavoring, as well as in traditional herbal formulas. It is part of the diet in China, Japan, and Asian communities in the United States. **Jarrow Formulas®** uses Red Yeast Rice that is carefully fermented to avoid the presence of *citrinin*, a naturally occurring mycotoxin. Coenzyme Q10 complements Red Yeast Rice in supporting healthy cardiovascular function.*

Keep out of the reach of children.

* These statements have not been evaluated by the Food and Drug Administration. This product is not intended to diagnose, treat, cure or prevent any disease.



Kaneka Q10® is a registered trademark of the Kaneka Corporation of Japan.

Lot #. Best Used Before:

Jarrow
FORMULAS

Red Yeast Rice

+ Co-Q10

Supports Cardiovascular Health*

120 DIETARY
VEGGIE CAPS SUPPLEMENT



www.nsfnon GMO.org

© 2019 Jarrow FORMULAS®

Supplement Facts

Serving Size 2 Capsules
Servings Per Container 60

	Amount Per Serving	% DV
Red Yeast Rice (<i>Monascus purpureus</i>)	1,200 mg	†
Coenzyme Q10 (Ubiquinone)	100 mg	†

† Daily Value not established.

Other Ingredients: Cellulose, magnesium stearate (vegetable source) and silicon dioxide. Capsule consists of hydroxypropylmethylcellulose.

No wheat, gluten, soybeans, dairy, egg, fish/shellfish, or peanuts/tree nuts.

Suitable for vegetarians/vegans.

Made in the U.S.A.

Store in a cool, dry place.

Distributed Exclusively by:
Jarrow FORMULAS®
Superior Nutrition and Formulation™
Los Angeles, CA 90035-4317
www.Jarrow.com
(866) 459-4154

06219RYR PROD # 133001



FILL LINE